



# C10CX

## Linear Color Doppler

[www.thaimedicaltrading.co.th](http://www.thaimedicaltrading.co.th)





Frequency

Depth

Mode

Applications

7.5-10MHz

20-100mm

B Mode  
B+M  
PW Doppler  
Color Doppler  
Power Doppler

MSK  
Breast  
Thyroid  
Superficial  
Vascular  
Pediatrics  
Carotid  
Small Parts



THAI MEDICALS TRADING Co.,Ltd.

67/17 MAHASAWAT PHUTTHAMONTHON NAKHONPATHOM

THAILAND 73170 TEL. 081-3451047,098-6329351

# Download our ultrasound software from App Store or Google Play



- ✓ Unlimited users
- ✓ Free software
- ✓ No hidden fees
- ✓ Support printing reports
- ✓ Chinese, English, Russian, Italian Spanish, Portuguese(Brazil)

## Software Controls

- GN(gain)
- D(depth)
- ENH(enhancement)
- DR(dynamic range)
- F(frequency)
- FocusPos
- PRF
- WF
- Mode
- 8TGC
- Biopsy
- Annote
- Page left/right up/down

## Measurement

- B: Length, Area/Circumference, Angle, Trace, Distance
- B+M: Heart Rate, Time, Distance
- B+PW: Velocity, Heart Rate(2), S/D, Depth



**THAI MEDICALS TRADING Co.,Ltd.**

67/17 MAHASAWAT PHUTTHAMONTHON NAKHONPATHOM

THAILAND 73170 TEL. 081-3451047,098-6329351



Package includes:

1\* wireless probe

1\* wireless charger

1\* USB charger

1\* wrist strap

1\* plug

1\* portable aluminium alloy box

## Connection Type

WiFi (802.11n/20MHz/5G)

Probe built-in WiFi signal

## Supporting Devices

Tablet, Smart phone, PC, etc

## Working system

Android, IOS, Windows

## Buttons

3 programmable buttons

## Lithium Battery

4200 mAh. 3.7V

Working time:  $\geq 3h$

Standby time:  $\geq 12h$

Battery life 5-8years

wireless charging

## Dimensions

158(h)x70(w)x30(d) mm

Net weight 223g

Elements 192

Footprint L45

## Boot-up

Platform dependent, generally <5 seconds




**THAI MEDICALS TRADING Co.,Ltd.**

67/17 MAHASAWAT PHUTTHAMONTHON NAKHONPATHOM

THAILAND 73170 TEL. 081-3451047,098-6329351

# Offer great images to everyone



THAI MEDICALS TRADING Co.,Ltd.

67/17 MAHASAWAT PHUTTHAMONTHON NAKHONPATHOM

THAILAND 73170 TEL. 081-3451047,098-6329351