



# C10RL


## Double head



# USB&Wi-Fi 3 in 1 probe


can be connected to a Windows/Android/IOS device  
(Tablet, smart phone, PC etc) via WiFi or USB cable.

- Primary Care
- Emergency
- EMS
- Veterinary
- Wellness
- Interventional
- Med School
- OB/GYN
- Anesthesia
- Breast
- Pain Management
- Sports Medicine
- Global Health
- Clinical Research



Freq:  
7.5/10MHz  
Depth:  
20-100mm  
Footprint:  
L25

**Linear**



Type-C

125mm  
N.W:90g

**Convex + Cardiac**


Freq:  
2.5/5MHz  
3.5/5MHz  
Depth:  
90-160mm  
90-305mm  
Radius:  
R45



You can scan anybody anytime anywhere conveniently, such as a clinic, hospital ward, emergency room, or operating theater.

Application: Abdomen、Obstetrics / Gynecology、Cardiac、Kidney、Urology、Breast、Thyroid、MSK、Vascular、Pediatrics、Carotid

Download our ultrasound software  
from App Store / Google Play



Probe built-in WiFi signal

Lithium Battery  
2800mAh

Working time:  $\geq 3h$   
Standby time:  $\geq 12h$

Modes: B, B+M,  
Color Doppler, PW,  
Power Doppler

Support working while charging

## Software Controls

- ✓ Support report edit and print
- ✓ Chinese, English, Russian, Italian, Spanish, Portuguese (Brazil)
- GN (gain), D (depth), ENH (enhancement), DR (dynamic range)
- F (frequency), Focus Pos, PRF, WF, Mode, 8TGC, Biopsy, Annotate,
- Page left/right up/down

## Measurement

Length, Area/Circumference, Angle, Trace, Distance GA (CRL, BPD, GS, FL, HC, AC)  
EFW (BPD, FL), Heart Rate, Time, Distance, Velocity, Heart Rate (2), S/D, Depth



THAI MEDICALS TRADING Co., Ltd.

67/17 MAHASAWAT PHUTTHAMONTHON NAKHONPATHOM  
THAILAND 73170 TEL. 081-3451047, 098-6329351


# Clinical images




Thyroid C Mode




Carotid C Mode



Vascular PW Mode



Liver B Mode



Cardiac C Mode



Kidney C Mode



**Anybody**

**Anywhere**

**Anytime**



**THAI MEDICALS TRADING Co.,Ltd.**

67/17 MAHASAWAT PHUTTHAMONTHON NAKHONPATHOM

THAILAND 73170 TEL. 081-3451047,098-6329351