

# EXLED300 Series



## Technical Parameter

Head dia	300mm
Illuminance	70,000-120,000Lux
Color temperature	4500±500K
Color reduction index(Ra)	>97
Size of light field	90-200mm
Service light of illuminant	50,000h
Power supply voltage	220V/50Hz
Rated power of bulb	3.2V/1w
Bulb	LED bulb
LED bulb quantity	24PCS



EXLED300/300



EXLED300



EXLED300EM



EXLED300M

## Product Functions Description

- 01/ Y-Shaped light head for easy positioning around physician's head



- 03/ Sterilize Handle



- 05/ Upper and down riser



- 07/ Stable base



- 02/ ON/OFF switch  
Illuminance adjustment  
Up to 120000 lux of pure, white bright light.



- 04/ 24 PCS Osram bulb, 50000h long service life



- 06/ Emergency panel  
Display the electricity and provide electrical socket.



- 08/ Medical silent castors with brakes

