

LED Operating Shadowless Lamps



CreLed 5500M Mobile LED Shadowless Operating Lamp



CreLed 3300M Mobile LED Shadowless Operating Lamp



Hollow CreLed 5500M Mobile LED Shadowless Operating Lamp

Type	CreLed 5500M	CreLed 3300M	Hollow CreLed 5500M
Lamp Head Diameter(mm)	500	600	500
Quantity of LED Bulbs	48	45	21
Illumination (Lux)	60,000-140,000	60,000-130,000	60,000-120,000
Color Rendering Index (Ra)	≥96	≥95	≥95
Color Temperature (K)	3500-5000	3500-5000	4000±500
Lighting Depth (mm)	≥1200	≥1200	≥1200
Spot Diameter (mm)	150-350	150-350	150-350
Dimming Range	12 levels	12 levels	12 levels
Operation Mode(Background, Surface, Deep Cavity and Normal Mode)	Yes	Yes	Yes
Bulb Service Life (hours)	≥50,000	≥50,000	≥50,000
Back-up Battery	≥8 hours	≥8 hours	≥8 hours



CreLed 150M LED Examination Lamp



(Mobile Type)



(Ceiling Type)

CreLed 310M LED Examination Lamp

Specifications	CreLed 150M	CreLed 310M
Illumination(Lux)	5,000-20,000	≥60,000
Color Temp	4000±500K	4000±500K
CRi(Ra)	≥93	≥93
Power Supply	110-130/220-240 V, 50/60 Hz	110-130/220-240 V, 50/60 Hz
LED Bulb Power(W)	3W/9V	42W/24V
Bulb Service Life	more than 50,000 hours	more than 50,000 hours
Bulb Number	1pc	14pcs
Highest Height of Head	190cm	280cm
Back-up Battery	N/A	4-6hours (Mobile Type)