

AQUARIUS

A Compact and Convenient
Patient Monitoring System



FEATURES



Mobility

Product dimension:
Unit: 173mm Wide x 152 mm Deep
Without Trolley: 274mm High
With Trolley: High Adjustable

Usability

Built-in rechargeable battery -
Standard 5 hour uninterrupted working
time, can be upgraded to last 15 hours!



Quality

8.4" TFT color LCD
Optional Touch Screen
High resolution 600*800.
Up to 6 parameter or 3 wave form display

Convenience

Multi-selectable user screen options

Introducing Aquarius

Aquarius from Northern Meditec is a compact, portable and simple to use vital signs monitor with an exciting array of accessories and options.

Get in touch with us today for more information...



Aquarius simplifies patient monitoring...

Optional Accessories

Touch screen, ECG, Quick Temp, Bar code scan, EtCO₂, Nellcor/Masimo SpO₂, Wired/Wireless CMS, Printer, Plus many more...



Touch Screen Interface



SpO₂ Sensors:
Neonatal, Paediatric &
Adult finger clip

Hard reusable sensor
Soft reusable sensor
SpO₂ Extension Cable



Infrared Tympanic
Thermometer Set



Reusable Cuffs:
Neonatal - 7-13cm
Infant - 12-19cm
Regular - 20 - 26cm
Adult - 25 - 35cm
Large Adult - 34 - 33cm



Height Adjustable Trolley
with Basket and Castors

FEATURES OF AQUARIUS

- Venipuncture mode, Spot check mode
- SpO₂ support PVI & PI, low perfusion 0.2%
- SpO₂ pulse-tone modulation (Pitch Tone)

Making more accurate screening and precise diagnoses

- Wired/Wireless CMS, support HL7
- Graphical & tabular trend review (120 hrs)
- Full disclosure waveforms review & storage

Clinical Applications:

Hemodynamic, Dose, Oxygenation, Ventilation, MEWS, Renal function



THAI MEDICALS TRADING CO.,LTD

Address : 67/17 Moo.1 Mahasawat, Phutthamonthon Nakhonpathom 73170 Thailand

Tel: +66 81-3451047 Web: www.tmt-med.com, Web: www.thaimedicaltrading.co.th E-mail : infinity_bio@hotmail.com