

# CM1200B


Electrocardiograph



## Features

- 5.6 inch LED screen
- Lead status indication
- Simple & clear keyboard layout, back light support
- Internal memory of 300 ECG records, built-in SD card (optional) can save 40000 cases
- High accurate digital filter (EMG filter, AC filter, Drift filter, lowpass filter)
- 120 seconds ECG waveforms review and print
- Support 210mm rolling, Z-folded paper
- USB port for external printer connection & data transmission
- PC-ECG management system
- Working mode: Auto, Manual, Rhythm

## Interface



## Various printing mode



Printing mode: 12X1, 6X2+1R, 3X4+1R



# CM1200B

Electrocardiograph



THAI MEDICALS TRADING CO.,LTD  
67/17 Moo.1 Mahasawat Phutthamonthon  
Nakhonpathom 73170 Tel: +66 81-3451047



# CM1200B

Electrocardiograph



Hidden handle




Alphanumeric keyboard;  
Back light for button



Memory of 300 ECG records;  
USB for data transmission;  
PC-ECG management system

## Clear professional report



12X1 + Average template

✓ 5.6 inch LED screen, easy for observation

✓ Internal aluminum alloy rack, shockproof, Anti-interference design

✓ High quality printing system, print is clear and smooth

✓ Internal modular design, easy for maintenance

