



Miller0



Mac3



Miller1



Mac4



Mac2

VIDEO LARYNGOSCOPE WITH REUSABLE BLADE MODEL NO.HYHJ-KC

Blade	Miller0	Miller1	Mac2	Mac3	Mac4
Application	Infant&S Child	M Child	S Adult&L Child	M adult	L Adult
Length(mm)	59.9	77.4	122	140	153
Height(mm)	18.55	18.35	29.23	29.23	29.23
Sterilization	*2.0% Glutaral fully soak 20 minutes. *Ethylene oxide sterilization 20 minutes *Low Temperature Plasma				





VIDEO LARYNGOSCOPE WITH DISPOSABLE BLADE MODEL:HYHJ-1320(NEW)

Blade	Mac1	Mac2	Mac3	Mac4	Mac5
Application	Pediatric	Child	M Adult	L adult	Difficult Airway
Length(mm)	106.14	116.14	126.14	132.14	141.14
Height(mm)	18.45	24.7	28.53	31.46	36.71
Interface	Interface:Type-C		Video Output:HDMI		

