



Hi-kV Portable X-ray **XUPERB**



Overview

- Clear Images (focal spot : 0.8mm)
- Light weight, perfect portable (3.8kg with Battery)
- Amazing performance (80kV / 5mA)
- More than 100 shots with full charge battery
- Able to shot X-ray even while charging
- Clear collimator beam
- Laser SID measurement




High Quality Images



Xuperb EZ-805 Specifications

X-ray Generator	Type	MEDICAL
	Power Output	400 W
X-ray Performance	Tube Voltage	50 ~ 80 kV (1 kV step)
	Tube Current	2 ~ 5 mA (1 mA step)
X-ray Control	Exposure Time Set	Time : 0.03 ~ 2.00 sec 0.03 ~ 0.20 (10 msec step) 0.20 ~ 2.00 (50 msec step)
Collimator	SID	SID 100cm / 43cm x 43cm
X-ray Tube	Anode Type	Stationary
	Anode Angle	16 Degrees
	Focal Spot	0.8 mm
	Anode Heat Capacity	8000 J
	Max. Anode Cooling Rate	110 W
	Nominal Anode Input Power at 0.1 s (DC)	1280 W
Power Requirement	Battery	Li-Polymer 22.2 VDC, 1450 mAh
	Number of Exposure	About more than 100 times at 70 kV, 5 mA, 0.1 sec
	AC Power Adapter	AC 100 ~ 240 V, 50/60 Hz
	Charging Time of Battery	About 1.5 ~ 2.0 hour
Weight	Main Body	3.8 kg
Dimension	Main Body	266 x 260 x 163 mm

Dimension



Various Option



Basic Stand



Deluxe Stand



Wall Bracket



Wall Stand



SIGNERS CO., LTD.

E-mail : supia@signers.co.kr

